



Trilogy EV300 *evolution of care*

simple, portable, reliable, adaptable

February 2020

innovation  you

Introduction


The Trilogy EV300 ventilator provides:

Invasive and non-invasive positive pressure ventilation for the care of patients ≥ 2.5 kg through adults.


The ventilator can measure, display, record, and alarm SpO_2 , FiO_2 , CO_2 , respiratory rate, and pulse rate data when integrated with the appropriate accessories.

The ventilator is suitable for use in the hospital transport setting*

*Not MRI compatible



Simple
Portable
Reliable
Adaptable



Simple

User-friendly platform


Patient-friendly performance

8" touchscreen



Simple

Quick access 100% oxygen flush button to deliver for 2 minutes




Simple

To prevent accidental therapy changes, use the **touchscreen lock**.

This is a temporary touchscreen lock, which can be changed back by tapping anywhere on the screen and following the onscreen instruction.

For automatic touchscreen lock, go to the Options screen then Device Options and select Automatic Touchscreen Lock On.



Simple

Quick start up. Ideal body weight (IBW) is calculated based on height and gender.

This information is used to establish default therapy and alarm settings, including tidal volume and alarms based on tidal volume. This information also limits setting ranges.



Simple

Adjust alarms as needed.



Simple

Modes and settings: **Trilogy 202 to Trilogy EV300**

Trilogy 202	Trilogy EV300	Description
AC	A/C - VC	Assist Control (Volume Control) mode provides volume-controlled mandatory or assist-control breaths. The set inspiratory time applies to all breaths.
CV		If you want to replicate CV mode where the ventilator triggers and cycles all breaths then set the trigger type to OFF.
PC	A/C - PC	Assist Control (Pressure Control) mode provides pressure-controlled mandatory or assist-control breaths. The set inspiratory time applies to all breaths. <i>Optional AVAPS.</i>
T		If you want to replicate T mode where the ventilator triggers and cycles all breaths then set the trigger type to OFF.
SIMV	SIMV-VC	Synchronized Intermittent Mandatory Ventilation (Volume Control) mode is similar to SIMV-PC, but with volume control.
PC-SIMV	SIMV-PC	Synchronized Intermittent Mandatory Ventilation (Pressure Control) mode is a pressure control mode that provides a mixture of mandatory, assist-control and spontaneous breaths with optional pressure support. It guarantees one mandatory breath in each cycle. The breath rate determines the length of the cycle. <i>Optional: Inspiratory Time min/max. for the spontaneous breaths.</i>

Simple

Modes and settings: **Trilogy 202 to Trilogy EV300**

Trilogy 202	Trilogy EV300	Description
S	PSV	Pressure Support Ventilation mode is patient-triggered, pressure-limited, and flow-cycled. The patient determines the breath rate and timing so it is recommended to set back-up ventilation. <i>Optional: AVAPS and Ti min/max.</i>
S/T	S/T	Spontaneous/Timed is a bi-level therapy mode where each breath is patient-triggered and patient-cycled, or ventilator-triggered and ventilator-cycled.
CPAP	CPAP	In Continuous Positive Airway Pressure mode , all breaths are spontaneous with the CPAP set pressure delivered in both inhalation and exhalation.
AVAPS-AE	AVAPS-AE	AVAPS-Auto EPAP mode automatically adjusts pressure support, to maintain the target tidal volume, and EPAP, to maintain a patent airway, within the set min/max ranges; and simplifies the set-up of the backup breath rate when set to auto. <i>Note: auto back-up rate maximum is 20bpm. Optional: Inspiratory Time min/max.</i>

Simple

Modes and settings: **Trilogy 202 to Trilogy EV300**

Trilogy 202	Trilogy EV300	Description
-	Inspiratory Time Min/Max	Once enabled, this setting treats inspiration time as a variable value for patient-initiated, patient-cycled breaths. It is available in S/T, PSV, SIMV-PC, SIMV-VC, and AVAPS-AE modes, under Advanced in the Prescription Settings window.
AVAPS Rate	AVAPS Speed	This sets the maximum rate of change in pressure between the min and max values while AVAPS is seeking a volume target.
-	PC Breath (AVAPS-AE)	Available in AVAPS-AE mode. When PC Breath is on, the set inspiratory time applies to all breaths.
Sigh	Sigh	In Trilogy EV300, available in A/C-VC mode under Advanced in the Prescription window. Sigh volume can be set between 1.5 – 2.5 times the set volume and the frequency every 50 – 250 breaths. While in Trilogy, sigh was fixed at 1.5 times the set volume every 100 breaths.
-	Back-up Ventilation	Available under Advanced in the Prescription window. When turned on an Apnea interval needs to be set in the alarm settings tab. Within the apnea interval; if no breaths are triggered by the patient, the vent delivers breaths at the set pressure of volume based on the Backup Rate and Backup Insp Time.

Simple Onscreen help

Entering a new prescription or placing a new circuit on the ventilator is simple thanks to the addition of onscreen help.


Simply tap the help icon (?) for more information regarding that prescription setting or alarm situation.



Simple Onscreen battery indicator

During ventilation you can check the battery status.

Tap the battery icons in the toolbar to see the status of each battery.



Simple
Portable
Reliable
Adaptable



Portable

Ultimate Portability

15 hours of battery.*

Hot swappable detachable battery provides uninterrupted therapy.**



*Nominal run time per method in International Electrotechnical Commission (7.5 hr/battery). Detachable battery charge time 0% to 80% is 2.5 hours, Internal battery charge time 0% to 100% is 3.5 hours. A/C-VC mode ActivePAP circuit, PEEP 3cmH₂O and Vt 800ml.

**When the internal battery is charged, batteries can be replaced without the ventilator pausing therapy.

Simple
Portable
Reliable
Adaptable



Reliable

Low Total Cost of Ownership



Trilogly EV300



Trilogly 202

Trilogly Evo Service Solution

Avg. 21 mins



Preventive maintenance time

Trilogly Service Solution

Avg. 1 hour 40 mins

1,200 cycles



Battery cycles

475 cycles

1 year*

*Inspection, cleaning and filter replacement only



Service interval

10,000 hours / 2 years

Simple
Portable
Reliable
Adaptable



Adaptable

Seamlessly transition across care environments utilizing the same clinical technology



Different care
settings

Same clinical
technology

Adaptable

Evolution of ventilator technology

- ✓ Oxygen blender, low flow oxygen connection and FiO_2 cell
- ✓ 5 prescriptions
- ✓ 4 circuits: single and dual limb
- ✓ Circuit Calibration
- ✓ Tubing Compliance Compensation
- ✓ T_i min/max
- ✓ Flow Trigger 0.5
- ✓ Rise Time 0
- ✓ Dynamic Parameters
- ✓ AVAPS updates
- ✓ AVAPS-AE updates
- ✓ End-tidal CO_2 (optional)

Adaptable Oxygen

FiO₂ sensor access
on back panel

Up to 30Lpm low flow O₂

Oxygen Blender Module



Adaptable Prescriptions

Program up to 5 Prescriptions (presets).

Example:

Patient on A/C-VC mode undergoes daily weaning trials on CPAP



Adaptable Circuits



Passive

Active PAP

Active Flow

Dual Limb

	Passive	Active PAP	Active Flow	Dual Limb
Infant (9-13mm)				✓
Ped (14-18mm)	✓	✓	✓	✓
Adult/Ped (19mm)	✓	✓	✓	✓
Adult (20-22mm)	✓	✓	✓	✓
Min Set VT	50 ml	50 ml	35 ml	35 ml
External Flow Sensor Required			✓	✓

Adaptable Tubing compliance compensation

Trilogy EV300 excludes any losses in tidal volume due to the circuit.

Trilogy EV300 includes a default calibration providing automatic tubing compensation for the recommended circuits in the accessory guide.



Adaptable Circuit calibration

Volume losses in circuit tubing can be calculated and programmed into the Trilogy EV300 using the calibration method.




Adaptable Ti min/max

Available in S/T, PSV,
SIMV-PC, SIMV-VC, and
AVAPS-AE modes

Access under **Advanced**

Applicable to spontaneous
breaths only

Control / mandatory breaths
retain set **Insp. Time**



Adaptable Flow trigger

Flow trigger can be set to 0.5 L/min to offer increased sensitivity for your weakest patients.



Adaptable Rise Time

Rise Time is now even faster than Trilogy, and can be set to 0 to adapt to the needs of your patients.

Note: You can tap on the Help icon whenever it is visible and a screen will appear for information concerning that section.



Adaptable

AVAPS

Available in A/C-PC, S/T, and PSV modes

AVAPS Speed

- Replaced AVAPS Rate (of change) on Trilogy

AVAPS Startup

- First minute not limited by Speed setting
- Next session starts with the previous sessions final inspiratory pressure




Adaptable

AVAPS

Available in A/C-PC, S/T, and PSV modes

Algorithm resets to pressure midpoint when:

- AVAPS restart icon  is tapped
- Changing to another pre-set prescription, then changing back

Algorithm does not reset to pressure midpoint when:

- Changing the target tidal volume
- Changing the insp. pressure ranges



Adaptable AVAPS-AE additional flexibility

PC Breath – On/Off



Adaptable AVAPS-AE additional flexibility

PS Min/Max can go to 0

Please note that PS Min/Max will change to PC Min/Max when PC Breath is set to On.



Adaptable AVAPS-AE additional flexibility


Automatic algorithm restart

- AVAPS restarts at pressure midpoint
- EPAP returns to EPAPmin for 100 breaths
- AutoBUR (if enabled) restarts




Adaptable

Dynamic parameters




Dyn C

Static Compliance of respiratory system (lungs + chest wall), measured dynamically. Ratio between the change in volume to the change in pressure.




Dyn R

Airway Resistance
Estimate of the change in pressure divided by the air flow through the airways.



Dyn P_{plat}

Plateau pressure is the maximum pressure applied to small airways and alveoli during positive-pressure mechanical ventilation.



AutoPEEP

Estimate of the pressure (above PEEP) that exists in the patient airway at the end of exhalation.

Adaptable Dynamic parameters



Passive

Active PAP

Active Flow

Dual Limb


	Passive	Active PAP	Active Flow	Dual Limb
Dynamic Parameters*	✓		✓	✓
A/C-VC	✓		✓	✓
A/C-PC	✓		✓	✓
SIMV-VC	✓		✓	✓
SIMV-PC	✓		✓	✓

* During ventilator-initiated mandatory breaths and patient-initiated mandatory breaths

Adaptable

Pediatric Trached Patient Example:

Pediatric patient with tracheostomy tube on Trilogy EV300 had an increase in resistance noted over a 300 second period that was resolved after suctioning.





Trilogy EV300



Trilogy 202

Intended Use (weight)	>2.5 kg patient intended use (15 mL pressure modes / 35 mL volume modes)	>5 kg patient intended use
Battery	~7.5 internal + ~7.5 detachable	~3 internal + ~3 detachable
Circuits	Passive, Active PAP, Active Flow, Dual Limb	Passive, Active PAP, Active Flow, (MPV)
Pre-sets	5 pre-set prescriptions	2 pre-set prescriptions
Standby	✓	✗
Modes	Pressure - CPAP, S/T, PSV, A/C-PC, SIMV-PC, AVAPS-AE Volume - A/C-VC, SIMV-VC	Pressure - CPAP, S, S/T, T, PC, PC-SIMV, AVAPS-AE, PC-MPV Volume - AC, CV, SMIV, AC-MPV
AVAPS	First minute not limited by speed setting	Always limited by rate of change setting
Set Pressure (max)	60 cmH ₂ O	30 cmH ₂ O
Ti Min/Max	Spont. breaths (S/T, PSV, SIMV-PC, SIMV-VC, and AVAPS-AE modes)	Only set Ti
Flow Trigger	0.5 – 9 Lpm	1 – 9 Lpm



Trilog EV300



Trilog 202

Rise Time	0 - 6	1 - 6
Backup Ventilation	✓	✗
Dynamic lung parameters with no insp/exp hold	Dyn C, Dyn R, P _{plat} , autoPEEP	✗
FiO₂ sensor	✓	✗
Enhanced monitoring	Waveforms, SpO ₂ , EtCO ₂	Waveforms
Memory/Data transfer	Internal Memory (2GB) Data Transfer via USB	No internal memory Data Transfer via SD card
Circuit compensation	Circuit and humidifier selection Circuit calibration (optional)	✗
Touch Screen GUI	Touch Screen GUI	Non-touch screen GUI
On screen Alarm Guidance	✓	✗
Service/Maintenance	1 year interval (Inspection, cleaning and filter replacement only)	10,000; 17,500; (alternating every 10K and 7.5K blower hrs)



Simple

Easy-to-learn user interface, configurable to the care environment



Portable

15 hours of battery life, mobile stand for easy transport, easily mounts on wheelchairs



Reliable

The most robust and durable device we've ever created



Adaptable

Stays with patients as their care settings and needs change

