



**PHILIPS**

Mother & Child Care

Fetal, maternal and  
newborn solutions catalog

# Your passion our commitment

Caring for mother and baby in labor and delivery,  
the NICU, and the well-baby nursery

[Go to overview >](#)



### Fetal monitoring supplies

Quality fetal monitoring supplies are essential to delivering accurate, reliable, in utero measurements.

3



### Maternal monitoring supplies

Choose from a broad portfolio of monitoring supplies proven for delivering accurate, reliable measurements during labor and delivery.

12



### Neonatal monitoring supplies

Philips offers a broad family of neonatal supplies designed specifically for monitoring fragile neonates in the NICU.

26



### Jaundice management

Proven jaundice management solutions based on our developmentally friendly approach to treating newborns.

53



### Positioning

Developmentally-appropriate positioning solutions that help support caregivers' efforts in caring for the most fragile infants.

57



### Calming and soothing

Innovative, evidence-based solutions designed to help comfort even the smallest babies.

66



### Therapy support products

Easy-to-use products that support a baby's therapy.

72



### Education and consulting programs

A broad portfolio of educational offerings dedicated to supporting your goals of developmental care for the most fragile infants.

78



## Fetal Monitoring supplies



> Fetal direct ECG monitoring	4
> Intrauterine pressure (IUP) monitoring	6
> Fetal belts	7
> Fetal paper	9
> Other fetal supplies	10



# Fetal direct ECG monitoring



## Philips fetal scalp electrodes

Product #	Description	Packaging	Use with:
1 989803137631 CE	Fetal spiral electrode, single spiral	25 electrodes	DECG legplate adapter cables: 989803137651, 989803137791, 989803137801, 989803137811;
2 989803137641 CE	Fetal spiral electrode, double spiral (available in Europe only)	25 electrodes	



# Fetal direct ECG monitoring



## Philips DECG legplate adapter cables

Product #	Description	Packaging	Use with:
1 CE	DECG reusable legplate adapter cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	Philips fetal monitors: M1350A, M1350B, M1350C, M1353A, M2703A, M2705A Transducers/modules: M1357A, M1364-60001, M2738A, M1365A, M2727A, M2735A, M2720A
2 CE	DECG reusable legplate adapter cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	GE Corometrics monitors: 115, 116, 118, 119, 173, 174; GE Series monitors: 120, 126, 128, 129, 250 using a GE Corometrics Y-adapter cable 1442AAO or socket adapter 1442BAO.
3 CE	DECG reusable legplate adapter cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	GE Corometrics monitors: 120, 126, 128, 129
4 CE	DECG reusable legplate adapter cable for use with 989803137631 and 989803137641 fetal spiral electrodes	1 cable	Spacelabs IM77 (other models of Spacelabs may require an additional adapter cable.)

## Attachment electrode

Product #	Description	Packaging	Use with:
5 CE	Attachment electrode, minimum shelf life 6 months	25 electrodes	DECG reusable legplate adapter cables: 989803137651, 989803137791, 989803137801, 989803137811



# Intrauterine pressure (IUP) monitoring (not available in USA)



## Koala™ IUP solution

Product #	Description	Packaging	Use with:
1 M1333A CE	Koala™ disposable IUP catheters	10 catheters	Fetal monitors: M1350A, M1350B, M1350C, M1353A, M2703A, M2705A
2 989803143931 CE	Reusable IUP connector cable (Avalon) for Koala™ disposable IUP catheter	1 cable	Patient module: M2738A; TOCO+ transducer: M2735A; Fetal monitor: M2703A, M2705A
3 M1334A CE	Reusable IUP connector cable for Koala™ disposable IUP catheter	1 cable	Fetal monitors: M1350A, M1350B, M1350C, M1353A




# Fetal belts




1



2



3



4

## Disposable belts

Product #	Description	Packaging	Use with:
1 M2208A €€	Disposable abdominal belt: yellow, with buttons, 1.3m (4'3") long, 60mm (2.36") wide	100 rolls	Philips and legacy fetal monitors, transducers, and modules: M1355A, M1356A, M1358A, M1364A, M1365A, M2720A, M2725A, M2726A, M2727A, M2734A, M2734B, M2735A, M2736A, M2738A, 866075, 866076, 866077

## Reusable belts

Product #	Description	Packaging	Use with:
2 M1562B €€	Reusable abdominal belt: gray, washable, water resistant, pre-cut, 1.3m (4'3") long, 50mm (1.96") wide	5 rolls	Philips and legacy fetal monitors, transducers, and modules: M1355A, M1356A, M1357A, M1358A, M1359A, M1364A, M1365A, M2720A, M2725A, M2726A, M2727A, M2734A, M2734B, M2735A, M2736A, M2738A
3 M4601A €€	Reusable leg belt: gray, washable, water resistant, bulk roll, 15m (49') long, 32mm (1.26") wide	1 bulk roll	Philips and legacy fetal monitors, transducers, and modules: M1357A, M1359A, M1364A, M1365A, M2738A
3 M4602A €€	Reusable abdominal belt: gray, washable, water resistant, pre-cut, 1.3m (4'3") long, 60mm (2.36") wide	5 rolls	Philips and legacy fetal monitors, transducers, and modules: M1355A, M1356A, M1358A, M1359A, M1364A, M1365A, M2720A, M2725A, M2726A, M2727A, M2734A, M2734B, M2735A, M2736A, M2738A, 866075, 866076, 866077
4 M4603A €€	Reusable abdominal belt: gray, washable, water resistant, bulk roll, 15m (49') long, 60mm (2.36") wide	1 bulk roll	



# Fetal belts




## Belt guides, buttons and transducer accessories

Product #	Description	Packaging	Use with:
1 989803143401	Belt clip for TOCO+ Plastic belt clip for TOCO+ transducer	6 clips	Avalon transducers: M2734A, M2734B, M2735A, M2736A
2 989803184851	Reusable Clips for use with Avalon CL transducers and fetal belts. (Clip allows the transducer to be positioned appropriately along the belt)	10 clips	Avalon cableless transducers: 866075, 866076, 866077
M1569A	Fetal Belt Buttons	10 buttons	Reusable belts: M1562B, M4601A, M4602A, M4603A
989803184841	Connector cap to cover DECG/IUP socket when not in use	10 connector caps	Avalon cableless transducers: 866075 and 866077
989803184861	Battery Replacement Kit	1x rechargeable lithium-ion battery, 4.55WH 2x Transducer X seal rings 1x Labeled tool lever 1x Tool body 1x Installation note	Avalon cableless transducers: 866075, 866076, 866077





# Fetal paper



## Chemical/thermal paper

Product #	Description	Packaging	Use with:
1 M1910A CE	Fetal recording paper, intended for use in the U.S.A.: orange grid, 1-min. (3cm) highlights, 0-100 mmHg IUP scale, 30-240 FHR scale.	40 packs (150 pages per pack)	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350A, M1350B, M1350C, M1351A, M1353A; Recorder: M2925A
2 M1911A CE	Fetal recording paper, intended for use in Europe: green grid, 0-100 mmHg IUP scale, 0-13 kPa scale, 50-210 FHR scale	40 packs (150 pages per pack)	
M1913A CE	Fetal recording paper, intended for use in Canada and Australia: green grid, 1-min. (3cm) highlights, 0-100 mmHg IUP scale, 50-210 FHR scale.	40 packs (150 pages per pack)	
M1913J	Fetal recording paper, intended for use in Japan: green grid, 1-min. (3cm) highlights, 0-100 mmHg IUP scale, 50-210 FHR scale	40 packs (150 pages per pack)	



## Other fetal supplies



1



2



3

### Miscellaneous fetal supplies

Product #	Description	Packaging	Use with:
1 989803143411 CE	Avalon FM30/50-remote event marker	1 event marker	Avalon fetal monitors: M2702A, M2703A, M2704A, M2705A
M4503A	Fetal record file kit: 5 sets of tabs, 60 record wrappers, 120 patient labels, 40 antepartum labels, 60 delivery labels, 1 file box label	1 kit	

### Modules, transducers, adapters

Product #	Description	Packaging	Use with:
2 M2738A CE	Avalon patient module: for IUP, fetal/maternal ECG	1 module	Avalon fetal monitors: M2702A, M2703A, M2704A, M2705A

### Ultrasound transmission gel

Product #	Description	Packaging	Use with:
3 40483A CE	Ultrasound transmission gel, 0.25L (8 1/2oz) refillable squeeze bottle Shelf life: at least 6 months	12 bottles	Philips and legacy ultrasound transducers: M1356A, M1358A, M2720A, M2726A, 866075, 866076, 866077
40483B CE	Ultrasound transmission gel, 5 liter (1.3gal) bottle use to refill 40483A. Shelf life: at least 6 months	1x 5 liter bottle	



# Other fetal supplies

## Fetal supplies kits

Product #	Description	Packaging	Use with:
862214	Standard antepartum kit, available in North America only	10x M1569A 40x M1910A 12x 40483A gel 100x M2208A belts 6x Toco+ clips Avalon Fetal	Avalon fetal monitors: M2702A, M2703A, M2704A, M2705A
862215	Standard antepartum kit, available in China, European Union, Latin America, Asia (except Japan) only	10x M1569A 40x M1911A 12x 40483A gel 5x M4602A belts 6x Toco+ clips	
862215	Standard antepartum kit, available in Japan	10x M1569A 40x M1913J 12x 40483A gel 5x M4602A belts 6x Toco+ clips	
862217	Intrapartum DECG kit (single spiral), not available in China, Mexico, Russia	25x 989803137631 FSE 25x 989803139771 electrodes	Avalon fetal monitors: M2703A, M2705A
862218	Intrapartum DECG kit (double spiral), available in the European Union only	25x 989803137641 FSE 25x 989803139771 electrodes	



## Maternal monitoring supplies



- > **Non-invasive blood pressure** 13
- > **Pulse oximetry** 20
- > **Maternal temperature** 24
- > **Miscellaneous supplies** 25



# Non-invasive blood pressure

## Reusable NIBP cuffs



### Comfort Care cuffs – single hose

Product #	Description	Packaging	Use with:
M1573A CE	Small adult, light blue tubing Bladder width: 10.5cm (4.1") Bladder length: 24cm (9.4") Limb circumference: 20.5cm to 28cm (8" to 11")	1 cuff	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1573XL CE	Small adult extra-long, light blue tubing with edge striping Bladder width: 10.5cm (4.1") Bladder length: 24cm (9.4") Limb circumference: 20.5cm to 28cm (8" to 11") Total cuff length: 52cm (20.5")	1 cuff	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1574A CE	Adult, navy blue tubing Bladder width: 13cm (5.1") Bladder length: 28cm (11") Limb circumference: 27cm to 35cm (10.6" to 13.8")	1 cuff	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1574XL CE	Adult extra-long, navy blue tubing with edge striping Bladder width: 13cm (5.1") Bladder length: 28cm (11") Limb circumference: 27cm to 35cm (10.6" to 13.8") Total cuff length: 66cm (26")	1 cuff	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1575A CE	Large adult, burgundy tubing Bladder width: 16cm (6.3") Bladder length: 35cm (13.8") Limb circumference: 34cm to 43cm (13.4" to 16.9")	1 cuff	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1575XL CE	Adult extra-long, navy blue tubing with edge striping Bladder width: 16cm (6.3") Bladder length: 35cm (13.8") Limb circumference: 34cm to 43cm (13.4" to 16.9") Total cuff length: 83.5cm (32.9")	1 cuff	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1576A CE	Thigh, gray tubing Bladder width: 20cm (7.9") Bladder length: 43cm (16.9") Limb circumference: 42cm to 54cm (16.5" to 21.3")	1 cuff	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



# Non-invasive blood pressure

## Reusable NIBP cuffs



### Easy Care cuffs – single hose

Product #	Description	Packaging	Use with:
M4554B/ M4554B5 CE	Small adult, light blue. Bladder width: 10.6cm (4.2") Bladder length 24cm (9.4") Limb circumference: 20.5cm to 28.5cm (8" to 11.2")	1 cuff/5 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4555B/ M4555B5 CE	Adult, navy blue. Bladder width: 13.5cm (5.3") Bladder length: 30cm (11.8") Limb circumference: 27.5cm to 36cm (10.8" to 14.1")	1 cuff/5 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4556B/ M4556B5 CE	Adult extra-long, navy blue. Bladder width: 13.5cm (5.3") Bladder length: 30cm (11.8") Limb circumference: 27.5cm to 36cm (10.8" to 14.1") Total cuff length: 66cm (26")	1 cuff/5 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4557B/ M4557B5 CE	Large adult, burgundy. Bladder width: 17cm (6.7") Bladder length: 38cm (15") Limb circumference: 35cm to 45cm (13.7" to 17.7")	1 cuff/5 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4558B/ M4558B5 CE	Large adult extra-long, burgundy. Bladder width: 17cm (6.7") Bladder length: 38cm (15") Limb circumference: 35cm to 45cm (13.7" to 17.7") Total cuff length: 52.5cm (20.7")	1 cuff/5 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4559B/ M4559B5 CE	Thigh, gray. Bladder width: 21cm (8.3") Bladder length: 48cm (18.9") Limb circumference: 44cm to 56cm (17.3" to 22")	1 cuff/5 cuffs	Fetal monitors: M1350B, M1350C, M2702A, M2703A, M2704A, M2705A, M8105AS, 863275, 863276, 863278



# Non-invasive blood pressure

## Reusable NIBP cuffs

### Reusable cuff assortment kits - Easy Care and Comfort Care

Product #	Description	Packaging	Use with:
864288	Easy Care adult kit – 4 sizes	Consists of 4 cuffs: M4554B, M4555B, M4557B, M4559B	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
864291	Easy Care adult long Kit – 6 sizes	Consists of 6 cuffs: M4554B, M4555B, M4556B, M4557B, M4558B, M4559B	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1578A	Comfort Care adult kit	Consists of 4 cuffs: M1573A, M1574A, M1575A, M1576A	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M1579XL	Comfort Care adult long kit	Consists of 3 cuffs: M1573XL, M1574XL, M1575XL	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



# Non-invasive blood pressure

## Reusable NIBP cuffs



### Limited-use, vinyl, multi-patient – Multi Care cuffs

Product #	Description	Packaging	Use with:
989803183331 CE	Small adult, light blue Bladder width: 10.5cm (4.1") Bladder length: 24cm (9.4") Limb circumference: 20.5cm to 28cm (8" to 11")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
989803183341 CE	Adult, navy blue Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27cm to 35cm (10.6" to 13.8")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
989803183351 CE	Adult extra-long, navy blue with edge striping Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27cm to 35cm (10.6" to 13.8") Total cuff length: 58cm (22.8")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
989803183361 CE	Large adult, burgundy Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference: 34cm to 43cm (13.4" to 16.9")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
989803183371 CE	Thigh, gray Bladder width: 7.9" (20cm) Bladder length: 16.9" (43cm) Limb circumference: 42cm to 54cm (16.5" to 21.3")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A





# Non-invasive blood pressure

## Single-patient NIBP cuffs



### Durable, soft material - Gentle Care cuffs

Product #	Description	Packaging	Use with:
M4574B CE	Small adult, light blue. Bladder width: 10.6cm (4.2") Bladder length: 24cm (9.4") Limb circumference: 20.5cm to 28.5cm (8" to 11.2")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4575B CE	Adult, navy blue. Bladder width: 13.5cm (5.3") Bladder length: 30cm (11.8") Limb circumference: 27.5cm to 36cm (10.7" to 14.1")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4576B CE	Adult extra-long, navy blue. Bladder width: 13.5cm (5.3") Bladder length: 30cm (11.8") Limb circumference: 27.5cm to 36cm (10.7" to 14.1") Total cuff length: 56cm (22")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4577B CE	Large adult, burgundy. Bladder width: 17cm (6.7") Bladder length: 38cm (15") Limb circumference: 35cm to 45cm (13.8" to 17.7")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4578B CE	Large adult extra-long, burgundy. Bladder width: 17cm (6.7") Bladder length: 38cm (15") Limb circumference: 35cm to 45cm (13.8" to 17.7") Total cuff length: 83.5cm (32.9")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
M4579B CE	Thigh, gray. Bladder width: 21cm (8.3") Bladder length: 48cm (18.9") Limb circumference: 44cm to 56cm (17.3" to 22")	10 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



# Non-invasive blood pressure

## Single-patient NIBP cuffs



### Limited durability, soft material - Single Care cuffs

Product #	Description	Packaging	Use with:
989803182291 CE	Small adult, light blue Bladder width: 10.6cm (4.2") Bladder length: 24cm (9.4") Limb circumference: 20.5cm to 28.5cm (8" to 11.2")	20 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
989803182301 CE	Adult, navy blue Bladder width: 13.5cm (5.3") Bladder length: 30cm (11.8") Limb circumference: 27.5cm to 36cm (10.7" to 14.1")	20 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
989803182311 CE	Adult extra-long, navy blue with edge striping Bladder width: 13.5cm (5.3") Bladder length: 30cm (11.8") Limb circumference: 27.5cm to 36cm (10.7" to 14.1") Total cuff length: 66cm (26")	20 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
989803182321 CE	Large adult, burgundy Bladder width: 17cm (6.7") Bladder length: 38cm (15") Limb circumference: 35cm to 45cm (13.8" to 17.7")	20 cuffs	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



# Non-invasive blood pressure

## NIBP air hoses



1



2

### Adult NIBP air hose

Product #	Description	Packaging	Use with:
1 M1598B CE	Adult pressure air hose, length: 1.5m (4.92')	1 air hose	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A
2 M1599B CE	Adult pressure air hose, length: 3.0m (9.84')	1 air hose	Fetal monitors: M2702A, M2703A, M2704A, M2705A, M1350B, M1350C, M1351A, M1353A



# Pulse oximetry

## Reusable SpO<sub>2</sub> sensors



1



2



3

### Reusable SpO<sub>2</sub> sensors

Product #	Description	Packaging	Use with:
1 M1191B CE	Adult finger glove sensor Length: 2m (6.5') Patient size: >50kg (>110lb) 8 pin, Philips connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A, (To add length, use optional SpO <sub>2</sub> extension cable: M1941A); Fetal monitors: M1350B, M1350C (Requires SpO <sub>2</sub> adapter cable: M1940A)
1 M1191BL CE	Adult finger glove sensor Length: 3m (9.8') Patient size: >50kg (>110lb) 8 pin, Philips connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A (Extension cable not allowed)
2 M1191T CE	Adult finger glove sensor Length: 45cm (1.5') Patient size: >50kg (>110lb) 9 pin, D-sub connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A, (Requires SpO <sub>2</sub> adapter cable: M1943A/AL); Fetal monitors: M1350B, M1350C (Requires SpO <sub>2</sub> adapter cable: M1900B)
3 M1192A CE	Pediatric/small-adult glove sensor Length: 1.5m (4.9') Patient size: 15 to 50kg (33 to 110lb) 8 pin, Philips connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A, (To add length, use optional SpO <sub>2</sub> extension cable: M1941A); Fetal monitors: M1350B, M1350C (Requires SpO <sub>2</sub> adapter cable: M1940A)
3 M1192T CE	Pediatric/small-adult glove sensor Length: 45cm (1.5') Patient size: 15 to 50kg (33 to 110lb) 9 pin, D-sub connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A, (Requires SpO <sub>2</sub> adapter cable: M1943A/AL); Fetal Monitors: M1350B, M1350C (Requires SpO <sub>2</sub> adapter cable: M1900B)



# Pulse oximetry

## Reusable SpO<sub>2</sub> sensors



### Reusable SpO<sub>2</sub> sensors (continued)

Product #	Description	Packaging	Use with:
1 M1194A CE	Pediatric/adult ear clip sensor Length: 1.5m (4.9') Patient size: >40kg (>88lb) 8 pin, Philips connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A, (To add length, use optional SpO <sub>2</sub> extension cable: M1941A); Fetal monitors: M1350B, M1350C (Requires SpO <sub>2</sub> adapter cable: M1940A)
2 M1196A CE	Pediatric/adult finger clip sensor Length: 3m (9.8') Patient size: >40kg (>88lb) 8 pin, Philips connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A (Extension cable not allowed)
M1196S CE	Adult finger clip sensor 2m (6,5') cable, Patient size: >40kg (>88lb), 8-pin, Philips connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A (To add length, use optional SpO <sub>2</sub> extension cable: M1941A)
3 M1196T CE	Pediatric/adult finger clip sensor Length: 90cm (2.9') Patient size: >40kg (>88lb) 9 pin, D-sub connector	1 sensor	Fetal monitors: M2703A, M2704A, M2705A, (Requires SpO <sub>2</sub> adapter cable: M1943A/AL); Fetal monitors: M1350B, M1350C (Requires SpO <sub>2</sub> adapter cable: M1900B)
M1627A CE	Wristband Velcro straps for SpO <sub>2</sub> sensor cables	10 straps	



# Pulse oximetry


## Single-patient SpO<sub>2</sub> sensors



1



2



3

### Single-patient SpO<sub>2</sub> sensors

Product #	Description	Packaging	Use with:
1 M1131A CE	Adult/pediatric finger clip sensor (Single-patient use) Length: 45cm (18") Patient size: >20kg (>44lb) 9-pin, D-sub connector	20 sensors	Fetal monitors: M2703A, M2704A, M2705A (Requires SpO <sub>2</sub> adapter cable: M1943A/AL); Fetal monitors: M1350B, M1350C (Requires SpO <sub>2</sub> adapter cable: M1900B) Adapter Cables: M1900B, M1943A, M1943AL SpO <sub>2</sub> Modules: M1020A and M1020B, Opt. A01 Multi-measurement servers: M3000A, M3001A, M3002A, Opt. A01 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863069, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283
2 M1133A CE	Neonatal/infant/adult Single-patient wrap sensor Length: 0.9m (35") Neonatal foot/hand, patient size: <3kg (<6.6lbs); Infant big toe/thumb, patient size: 10 to 20kg (22 to 44lbs); Adult finger except thumb, patient size: > 40kg (> 88lbs)	20 sensors	
3 M1134A CE	Adhesive-free neonatal/infant/adult Single-patient wrap sensor Length: 0.9m (35") Neonatal foot/hand, patient size: <3kg (<6.6lbs); Infant big toe/thumb, patient size: 10 to 20kg (22 to 44lbs); Adult finger except thumb, patient size: > 40kg (> 88lbs)	20 sensors	



# Pulse oximetry

## SpO<sub>2</sub> cables



1



2

### SpO<sub>2</sub> cables

Product #	Description	Packaging	Use with:
M1900B CE	SpO <sub>2</sub> adapter cable Length: 3m (9.8') Adapts 9-pin sensors to 12-pin sockets	1 cable	Fetal monitors: M1350B, M1350C
M1940A CE	SpO <sub>2</sub> adapter cable Length: 2.0 m (6.6') Adapts 8-pin sensors to 12-pin sockets	1 cable	Fetal monitors: M1350B, M1350C
1 M1941A CE	SpO <sub>2</sub> extension cable, Length: 2m (6.6') For 8-pin sensors to 8-pin, Philips sockets Adds additional length to sensor cable	1 cable	Fetal monitors: M2703A, M2704A, M2705A
2 M1943A CE	SpO <sub>2</sub> adapter cable (8-pin) Length: 1.1m (3.3') Adapts 9-pin sensors to 8-pin, Philips sockets	1 cable	Fetal monitors: M2703A, M2704A, M2705A
2 M1943AL CE	SpO <sub>2</sub> adapter cable (8-pin) Length: 3m (9.8') Adapts 9-pin sensors to 8-pin, Philips sockets	1 cable	Fetal monitors: M2703A, M2704A, M2705A



# Maternal temperature



## Maternal temperature

Product #	Description	Packaging	Use with:
1 989803179381	Genius 2 probe cover (non-CE marked)	6 bars x 16 covers each in one box, 22 boxes in 6 x 16 x 22 = 2112 probes covers.	Fetal monitors: M2702A, M2703A, M2704A, M2705A
1 989803179611 CE	Genius 2 probe cover (CE marked)	6 bars x 16 covers each in one box, 22 boxes in 6 x 16 x 22 = 2112 probes covers.	
2 989803180831	Genius 2 tympanic temperature probe	1 per box	





# Miscellaneous supplies



1



2

## Disposable ECG electrodes

Product #	Description	Packaging	Use with:
1 40493D CE	Adult foam ECG electrode	5 per pack 300 per case	Fetal/maternal monitors: M1350B, M1350C, M2702A, M2703A, M2704A, M2705A, M2720A
1 40493E CE	Adult foam ECG electrode	30 per pack (300 per case)	Fetal/maternal monitors: M1350B, M1350C, M2702A, M2703A, M2704A, M2705A, M2720A

## ECG Lead Set

Product #	Description	Packaging	Use with:
2 M1363A CE	MECG cable with reusable leadset for maternal ECG	1 cable	Fetal monitors/patient modules: M1364A, M1365A, M2727A, M2738A, M2735A, M2702A, M2703A, M2704A, M2705A, M2720A



## Neonatal monitoring supplies



> Capnography	27
> Spirometry	30
> Transcutaneous gas monitoring	31
> ECG cables and leads	33
> NeoLead ECG electrodes	35
> ECG electrodes	37
> Non-invasive blood pressure	40
> Pulse oximetry	47
> Temperature	51



# Capnography

## Sidestream CO<sub>2</sub>



### Microstream® FilterLine®

Product #	Description	Packaging	Use with:
1 M1923A CE	FilterLine H set, infant/neonatal, intubated, long term.	25 sets	Microstream ext. server: M3015A, M3015B; Monitors: 863266, 863066, M8105A, M8105AS
2 989803160261 CE	FilterLine® H set, infant/neonatal, long	25 sets	Microstream ext. server: M3015A, M3015B; Monitors: 863266, 863066, M8105A, M8105AS
3 989803159581 CE	VitaLine H set, infant/neonatal, intubated, long term, high ambient humidity environments.	25 sets	Microstream ext. server: M3015A, M3015B; Monitors: 863266, 863066, M8105A, M8105AS
4 M4691A CE	CapnoLine H, infant/neonatal	25 sets	Microstream ext. server: M3015A, M3015B; Monitors: 863266, 863066, M8105A, M8105AS
989803178001 CE	CapnoLine H O <sub>2</sub> , infant/neonatal	25 sets	Microstream ext. server: M3015A, M3015B; Monitors: 863266, 863066, M8105A, M8105AS
989803178011 CE	CapnoLine, infant/neonatal, long	25 sets	Microstream ext. server: M3015A, M3015B; Monitors: 863266, 863066, M8105A, M8105AS
989803178021 CE	CapnoLine, infant/neonatal	25 sets	Microstream ext. server: M3015A, M3015B; Monitors: 863266, 863066, M8105A, M8105AS



# Capnography

## Sidestream CO<sub>2</sub>




### LoFlo sidestream CO<sub>2</sub>

	Product #	Description	Packaging	Use with:
1	M2741A CE	Sidestream CO <sub>2</sub> sensor	1 sensor	Capnography extension server : M3014A; Monitors: M3002A, M8102A, M8105A
2	M2746A CE	CO <sub>2</sub> nasal cannula, infant	10 sets	Sidestream sensor with capnography ext. server: M3014A; Monitors: M8105A, M8102A, M3002A
3	989803144471 CE	CO <sub>2</sub> /O <sub>2</sub> nasal cannula, infant/neonatal	10 sets	Sidestream sensor (M2741A) with capnography ext. server: M3014A; Monitors: M8105A, M8102A, M3002A
4	M2776A CE	Straight sample line, non-humidified	10 sets	Sidestream sensor (M2741A) with capnography ext. server: M3014A; Monitors: M8105A, M8102A, M3002A
5	M2777A CE	Straight sample line H, humidified	10 sets	Sidestream sensor (M2741A) with capnography ext. server: M3014A; Monitors: M8105A, M8102A, M3002A



# Capnography

## Mainstream CO<sub>2</sub>



### Mainstream CO<sub>2</sub> monitoring

Product #	Description	Packaging	Use with:
1 M2501A CE	Mainstream CO <sub>2</sub> sensor	1 sensor	Capnography extension server: M3014A; Monitors: M8105A, M8102A, M3002A,
2 M2516A CE	Airway adapter, reusable, infant ET tube size: less than 4mm, deadspace less than 1cc	1 airway adapter	Mainstream sensor (M2501A) with capnography ext. server: M3014A; Monitors: M8105A, M8102A, M3002A
3 M2536A CE	Airway adapter Single-patient use, infant ET tube size: less than 4mm, deadspace less than 1cc	10 airway adapters	M3014A sensor (M2501A) with capnography ext. server: M3014A; Monitors: M8105A, M8102A, M3002A



# Spirometry



## Flow sensors and combined CO<sub>2</sub>/flow

Product #	Description	Packaging	Use with:
1 M2782A CE	Combined CO <sub>2</sub> /flow sensor, infant/ neonatal, disposable ET: less than 4mm	10 sensors	Mainstream sensor (M2501A) with spirometry module (M1014A)
2 M2786A CE	Flow sensor, infant/ neonatal, disposable ET: less than 4mm	10 sensors	Spirometry module: M1014A



# Transcutaneous gas monitoring



## tcpO<sub>2</sub>/tcpCO<sub>2</sub>

Product #	Description	Packaging	Use with:
M1918A €€	Solid state tcpO <sub>2</sub> /tcpCO <sub>2</sub> transducer	Includes: 1 transducer, 1 cap	Module: M1018A
1 M1918B €€	Solid state tcpO <sub>2</sub> /tcpCO <sub>2</sub> transducer	Includes: 1 transducer, 1 cap, 1 calibration ring	Module: M1018A
2 989803151951 €€	tcpO <sub>2</sub> /tcpCO <sub>2</sub> membrane kit	10 membranes, 10ml electrolyte, 2 sensor caps, 5 cleaning papers	Transducer: M1918B
3 989803151961 €€	tcpO <sub>2</sub> /tcpCO <sub>2</sub> fixation kit	100 fixation rings, 4 bottles contact fluid @ 20ml	Transducer: M1918B
4 989803160371	Calibration ring	1 ring	Module: M1018A
15209-60010 €€	Accessory kit for tcpO <sub>2</sub> /tcpCO <sub>2</sub> transducer	10ml electrolyte, 12 membranes, 2 transducer caps, 25 cleaning papers, Instructions, O-ring remover	Transducer: M1918A
15209-60020 €€	Fixation kit for tcpO <sub>2</sub> /tcpCO <sub>2</sub> transducers	100 fixation rings, 4 (20ml) bottles contact liquid, instructions	Transducer: M1918A



# Transcutaneous gas monitoring

## tcpO<sub>2</sub>/tcpCO<sub>2</sub> (continued)

Product #	Description	Packaging	Use with:
15210-60010	Cal 1 gas cylinders for calibration (not available in Europe or Japan)	6 cylinders	Module: M1018A Calibration unit: 15210B Transducer: M1918A, M1918B Microstream extension server: M3015A, M3015B
15210-60020	Cal 1 gas cylinders for calibration (not available in Europe or Japan)	6 cylinders	
15210-64010	Gas Cal 1 gas cylinders for tcpO <sub>2</sub> /tcpCO <sub>2</sub> calibration (meets international shipping requirements)	6 cylinders	
15210-64020	Gas Cal 2 gas cylinders for tcpO <sub>2</sub> /tcpCO <sub>2</sub> calibration (meets international shipping requirements)	6 cylinders	
M2205A	Calibration tubes used to connect 15210B calibration unit to M1918A or M1918B transducer during calibration	5 tubes	Module: M1018A Calibration unit: 15210B Transducer: M1918A, M1918B





# ECG cables and leads



## 3-lead IntelliVue ECG patient cable system

- Lightweight and comfortable for tiny patients

Product #	Description	Packaging	Use with:
1 M1669A CE	3-lead ECG trunk cable, AAMI/IEC, 12-pin connector Length: 2.7m (9')	1 cable	Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 863088, 863317; Efficia Monitors 863300, 863301, 863302, 863303, 863304, 863320, 863321, 863322; Defibrillators M3534A, M3535A, M3536A (Heart Start MRx), M4735A (HeartStart XL), 861290 HeartStart XL+, DFM 100
2 989803170171 CE	3-lead OR ECG trunk cable, AAMI/IEC, 12-pin connector Length: 2.7m (9') (orange)	1 cable	Multi-Measurement Modules: M3001A, M3002A; MP2: M8102A; MP5: M8105A; SureSigns Monitors: 863064, 863065, 863066, 863068, 863317; Efficia Monitors 863300, 863301, 863302, 863303, 863304; Defibrillator: DFM 100
3 M1622A	3-lead miniclip lead set, AAMI Length: 0.5m (1.5'), grabbers, non-shielded (purple)	1 lead set	Trunk cable: M1669A and 989803170171
4 M1624A	3-lead miniclip lead set, AAMI Length: 0.7m (2.3'), grabbers, non-shielded (purple)	1 lead set	Trunk cable: M1669A and 989803170171
5 M1626A CE	3-lead miniclip lead set, IEC Length: 0.7m (2.3'), grabbers, non-shielded (purple)	1 lead set	Trunk cable: M1669A and 989803170171



# ECG cables and leads




## 5-lead IntelliVue ECG patient cable system

- Lightweight and comfortable for tiny patients

Product #	Description	Packaging	Use with:
1 M1668A CE	5-lead ECG trunk cable, AAMI/IEC, 12-pin connector Length: 2.7m (9')	1 cable	Modules: M1001B, M1002B; Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A; SureSigns Monitors: 863064, 863065, 863066, 863068, 863086, 863087, 863088, 863317; Efficia Monitors 863300, 863301, 863302, 863303, 863304, 863320, 863321, 863322; Defibrillators: M3534A, M3535A, M3536A (Heart Start MRx), M4735A (HeartStart XL), 861290 HeartStartXL+, DFM 100; Xper Flex Cardio: 860338
989803170181 CE	5-lead ECG OR trunk cable, AAMI/IEC, 12-pin connector Length: 2.7m (9') (orange)	1 cable	Multi-Measurement Modules: M3001A, M3002A; MP2: M8102A; MP5: M8105A; SureSigns Monitors: 863064, 863065, 863066, 863068, 863317; Efficia Monitors: 863300, 863301, 863302, 863303, 863304, 863320, 863321, 863322; Defibrillator: DFM 100
2 M1647A	5-lead Miniclip lead set, AAMI Length: 0.7m/1.3m (2.3'/4.2'), grabbers, non-shielded (purple)	1 lead set	Trunk cable: M1668A and 989803170181
3 M1648A CE	5-lead Miniclip lead set, IEC Length: 0.7m/1.3m (2.3'/4.2'), grabbers, non-shielded (purple)	1 lead set	Trunk cable: M1668A and 989803170181



# NeoLead ECG electrodes




## NeoLead electrodes (radiolucent)

- Lightweight and comfortable for tiny patients
- Long-lasting

Product #	Description	Packaging	Use with:
1 989803183081	Mini NeoLead, metallic pre-attached leadwire (24") pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. AAMI, Standard DIN connector. 3/4" round.	3 electrodes per pouch, 50 per box	ECG modules: M1001B, M1002B; Multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863064, 863065, 863066, 863068, 860338, 860335 Requires ECG trunk cable: M1668A or 989803170181
2 989803183091 CE	Mini NeoLead, metallic pre-attached leadwire (60cm) pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. IEC, Standard DIN connector. 19mm round.	3 electrodes per pouch, 50 per box	
3 989803183101	NeoLead, metallic pre-attached leadwire (24") pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. AAMI, Standard DIN connector. 1" round.	3 electrodes per pouch, 50 per box	
4 989803183111 CE	NeoLead, metallic pre-attached leadwire (60cm) pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. IEC, Standard DIN connector. 25mm round.	3 electrodes per pouch, 50 per box	



# NeoLead ECG electrodes




## NeoLead electrodes

- Lightweight and comfortable for tiny patients

Product #	Description	Packaging	Use with:
1 989803183121	NeoLead, radiolucent pre-attached leadwire (24") pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. AAMI, Standard DIN connector. 1" round.	3 electrodes per pouch, 50 per box	ECG modules: M1001B, M1002B; Multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863064, 863065, 863066, 863068, 860338, 860335 Requires ECG trunk cable: M1668A or 989803170181
2 989803183131 CE	NeoLead, radiolucent pre-attached leadwire (60cm) pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. IEC, Standard DIN connector. 25mm round.	3 electrodes per pouch, 50 per box	
3 989803183141	Micro NeoLead, radiolucent pre-attached leadwire (24") pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. AAMI, Standard DIN connector. 1" x 3/8".	3 electrodes per pouch, 50 per box	
4 989803183151 CE	Micro NeoLead, radiolucent pre-attached leadwire (60cm) pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. IEC, Standard DIN connector. IEC, 25mm x 9.5mm	3 electrodes per pouch, 50 per box	



# ECG electrodes




## Metallic, disposable lead wire electrodes

Product #	Description	Packaging	Use with:
1 13952A	Preattached metallic 46cm (18") leadwire electrode, non-shielded, AAMI solid gel. Dimension: 22mm x 33mm (7/8" x 1-1/4") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	Trunk cables M1669A and 989803170171
13952D CE	Preattached metallic 46cm (18") leadwire electrode, non-shielded, IEC solid gel. Dimensions: 22mm x 33mm, (7/8" x 1-1/4") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	
13953A	Preattached metallic 45cm (18") leadwire electrode, non-shielded, AAMI solid gel. Dimension: 22mm x 22mm (7/8" x 7/8") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	
1 13953D CE	Preattached metallic 46cm (18") leadwire electrode, non-shielded, IEC solid gel. Dimensions: 22mm x 22mm, (7/8" x 7/8") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	



# ECG electrodes



1


## Radiolucent, disposable lead wire electrodes

Product #	Description	Packaging	Use with:
1 13952B	Preattached radiolucent 46cm (18") leadwire electrode, non-shielded, AAMI solid gel. Dimension: 22mm x 33mm (7/8" x 1-1/4") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	Trunk cables M1669A and 989803170171
13952E CE	Preattached radiolucent 46cm (18") leadwire electrode, non-shielded, IEC solid gel. Dimensions: 22mm x 33mm, (7/8" x 1-1/4") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	
13953B	Preattached radiolucent 46cm (18") leadwire electrode, non-shielded, AAMI solid gel. Dimension: 22mm x 22mm (7/8" x 7/8") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	
13953E CE	Preattached radiolucent 46cm (18") leadwire electrode, non-shielded, IEC solid gel. Dimensions: 22mm x 22mm, (7/8" x 7/8") (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	



# ECG electrodes

Disposable snap connector electrodes



1



2



3

Product #	Description	Packaging	Use with:
1 989803192541 CE	Disposable, soft cloth backing, solid gel electrode. snap, grabber, 35mm (1.4") round. Shelf life: 24 months	30 electrodes/pouch 10 pouches/case	Pediatric/neonatal monitors
2 13951C CE	Disposable, repositionable solid gel electrode, snap, shelf life: 6 months. (3 electrodes mounted on plastic card) Dimensions: 22mm x 33mm (7/8" x 1-1/4")	30 electrodes/pack 300/case	Pediatric/neonatal monitors
3 13955C CE	Disposable, repositionable solid gel electrode, snap, shelf life: 6 months. (3 electrodes mounted on plastic card) Dimensions: 22mm x 22mm (7/8" x 7/8")	30 electrodes/pack 300/case	Pediatric/neonatal monitors



# Non-invasive blood pressure

## Reusable NIBP cuffs



1



2

### Reusable NIBP cuffs

Product #	Description	Packaging	Use with:
1 M1571A CE	Multi-patient Comfort Care cuff, 1 hose, infant Bladder width: 5.5cm (2.2") Bladder length: 15cm (5.9") Limb circumference: 10cm to 15cm (3.9" to 5.9")	1 cuff	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 866077, M8105AS, 860335, 860338, 863275, 863276, 863278
2 M4552B CE	Easy Care cuff, 1 hose, infant Bladder width: 5.5cm (2.2") Bladder length: 15cm (5.9") Limb circumference: 10cm to 15cm (3.9" to 5.9") Cuff length: 22cm (8.7")	1 cuff	NIBP module: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 866077, M8105AS, 863275, 863276, 863278
2 M4552B5 CE	Easy Care cuff kit, 1 hose, infant Bladder width: 5.5cm (2.2") Bladder length: 15cm (5.9") Limb circumference: 10cm to 15cm (3.9" to 5.9") Cuff length: 22cm (8.7")	5 cuffs	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 866077, M8105AS





# Non-invasive blood pressure

## Single-patient NIBP cuffs



1

### Soft Single-patient neonatal/infant cuffs

Product #	Description	Packaging	Use with:
1 M1866S CE	Neonatal, soft, Single-patient cuff, size #1, single hose	40 cuffs per box	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
1 M1868S CE	Neonatal soft Single-patient cuff, size #2, single hose	40 cuffs per box	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
1 M1870S CE	Neonatal soft Single-patient cuff, size #3, single hose	40 cuffs per box	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
1 M1872S CE	Neonatal soft Single-patient cuff, size #4, single hose	40 cuffs per box	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278



# Non-invasive blood pressure

## Single-patient NIBP cuffs

### Soft Single-patient neonatal/infant cuffs (continued)

Product #	Description	Packaging	Use with:
1 M1873S CE	Infant, soft, Single-patient cuff, size #5, single hose	40 cuffs per box	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS, 863275, 863276, 863278
1 989803167551 CE	Neonatal, soft, Single-patient, cuffs, single hose (consists of 5 x M1866S, 10 x M1868S, 20 x M1870S, 10 x M1872S & 5 x M1873S)	50 cuffs per box	NIBP modules: M1008B; Multi-measurement server: M3000A, M3001A, M3002A; Patient monitor: M8105A, M8102A, 863055, 863056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863066, 863068, 863069, 836070, 836071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, 836086, 863087, 863088, M8105AS



# Non-invasive blood pressure

## Single-patient NIBP cuffs



1

### Single-patient, vinyl, neonatal/infant cuffs

Product #	Description	Packaging	Use with:
1 M1866B CE	Neonatal Single-patient cuff, size 1 Circumference: 3.1cm to 5.7cm (1.2" to 2.2"), 1 hose. New connector configuration, replaces the M1866A.	40 cuffs	NIBP module: M1008B; Multi-Measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082; Legacy HP monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
1 M1868B CE	Neonatal Single-patient cuff, size 2 Circumference: 4.3cm to 8.0cm (1.7" to 3.1"), 1 hose. New connector configuration, replaces the M1868A.	40 cuffs	NIBP module: M1008B; Multi-Measurement servers: M3000A, M3001A, M3002A; Patient monitors M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082; Legacy HP monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
1 M1870B CE	Neonatal Single-patient cuff, size 3 Circumference: 5.8cm to 10.9cm (2.3" to 4.3"), 1 hose. New connector configuration, replaces the M1870A.	40 cuffs	NIBP module: M1008B; Multi-Measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082; Legacy HP monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
1 M1872B CE	Neonatal Single-patient cuff, size 4 Circumference: 7.1cm to 13.1 cm (2.8" to 5.2"), 1 hose. New connector configuration, replaces the M1872A.	40 cuffs	NIBP module: M1008B; Multi-Measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082; Legacy HP Monitors: 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278



# Non-invasive blood pressure

## Single-patient NIBP cuffs



2

### Single-patient neonatal/infant cuffs (continued)

Product #	Description	Packaging	Use with:
1 M1873B CE	Infant, single patient cuff, size 5 Circumference: 10.0cm to 15.0cm, (3.9" to 5.9") 1 hose New connector configuration.	40 cuffs	NIBP module: M1008B; Multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors; 78352A, 78352C, 78834C, M8105AS, 863275, 863276, 863278
1 989803167541 CE	Neonatal Single-patient NIBP kit contains 50 cuffs. Incorporates new connector configuration and replaces M1820-60020.	5x M1866B size 1, 10x M1868B size 2, 20x M1870B size 3, 10x M1872B size 4, 5x M1873B size 5	NIBP module: M1008B; Multi-measurement servers: M3000A, M3001A, M3002A, Patient monitors: M8105A, M8102A, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, Legacy HP Monitors; 78352A, 78352C, 78834C, M8105AS

### Gentle Care Single-patient NIBP cuffs

Product #	Description	Packaging	Use with:
2 M4572B CE	Soft Single-patient NIBP cuff, infant	10 cuffs	NIBP modules: M1008B; Multi-measurement servers: M3000A, M3001A, M3002A Patient monitors: M8105A, M8102A, 862474, 862478, 863055, 853056, 863057, 863058, 863059, 863060, 863061, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M3535A, M3536A, M8102A, M8105A, M8105AS, M8105AT, 860335, 860338, 863275, 863276, 863278



# Non-invasive blood pressure

## NIBP hoses



### Blood pressure air hoses

Product #	Description	Packaging	Use with:
1 M1596C CE	Neonatal blood pressure air hose. Length: 1.5m (4.92') Connects to Philips neonatal Single-patient cuffs only and incorporates new connector configuration. This cable is used for close patient proximity and replaces the M1596B.	1 air hose	NIBP modules: M1008B, Multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278
2 M1597C CE	Neonatal blood pressure air hose. Length: 3.0m (9.84') Connects to Philips neonatal Single-patient cuffs only and incorporates new connector configuration. This cable is used for classic bedside monitoring and replaces the M1597B.	1 air hose	NIBP modules: M1008B, Multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863079, 863080, 863081, 863082, M8105AS, 863275, 863276, 863278



# Non-invasive blood pressure

## NIBP hoses



1



2

### Blood pressure connectors

Product #	Description	Packaging	Use with:
1 989803167511	Connector upgrade kit for M1597B. Upgrade existing NIBP airhoses with an inner diameter of 2mm. Contains 20 female connectors with 2mm barb and instructions.	20 connectors	Upgrades M1597B to M1597C
2 989803167521	Connector upgrade kit for M1596B. Upgrade existing NIBP airhoses with an inner diameter of 3mm. Contains 20 female connectors with 3mm barb and instructions.	20 connectors	Upgrades M1596B to M1596C



# Pulse oximetry

## Single-patient SpO<sub>2</sub> Sensors



### Single-patient SpO<sub>2</sub> sensors

Product #	Description	Packaging	Use with:
1 M1132A CE	Infant Single-patient wrap sensor infant finger, Length: 0.9m (35") Patient size: >3 to 10kg (>6.6 to 22lbs.)	20 sensors	Adapter Cables: M1900B, M1943A, M1943AL SpO <sub>2</sub> Modules: M1020A and M1020B, Opt. A01 Multi-measurement servers: M3000A, M3001A, M3002A, Opt. A01 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863069, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283
2 M1133A CE	Neonatal/infant/adult Single-patient wrap sensor Length: 0.9m (35") Neonatal foot/hand, patient size: <3kg (<6.6lbs); Infant big toe/thumb, patient size: 10 to 20kg (22 to 44lbs); Adult finger except thumb, patient size: > 40kg (> 88lbs)	20 sensors	
3 M1134A CE	Adhesive-free neonatal/infant/adult Single-patient wrap sensor Length: 0.9m (35") Neonatal foot/hand, patient size: <3kg (<6.6lbs); Infant big toe/thumb, patient size: 10 to 20kg (22 to 44lbs); Adult finger except thumb, patient size: > 40kg (> 88lbs)	20 sensors	



# Pulse oximetry

## Reusable SpO<sub>2</sub> sensors



### Reusable SpO<sub>2</sub> sensors

Product #	Description	Packaging	Use with:
1 M1193A CE	Reusable neonatal wrap sensor Length: 1.5m (4.9') Neonatal foot/hand, patient size: 1kg to 4kg (2.2lbs to 8.8lbs) 8-pin connector	1 sensor	Adapter and Extension Cables: M1940A and M1941A SpO <sub>2</sub> Modules: M1020A and M1020B, Opt. A01 Multi-measurement servers: M3000A, M3001A, M3002A, Opt. A01 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863069, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283
2 M1193T CE	Reusable neonatal wrap sensor Length: 0.9m (35") Neonatal foot/hand, patient size: 1kg to 4kg (2.2lbs to 8.8lbs) 9-pin (D-sub) connector.	1 sensor	Adapter Cables: M1900B, M1943A, M1943AL SpO <sub>2</sub> Modules: M1020A and M1020B, Opt. A01 Multi-measurement servers: M3000A, M3001A, M3002A, Opt. A01 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863069, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283
3 M1195A CE	Reusable infant finger glove sensor Length: 1.5m (4.9') Infant finger, patient size: 4kg to 15kg (8.8lbs to 33lbs) 8-pin connector	1 sensor	Adapter and Extension Cables: M1940A and M1941A SpO <sub>2</sub> Modules: M1020A and M1020B, Opt. A01 Multi-measurement servers: M3000A, M3001A, M3002A, Opt. A01 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863069, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283





# Pulse oximetry

## SpO<sub>2</sub> adapter cables



### SpO<sub>2</sub> adapter cables

Product #	Description	Packaging	Use with:
1 M1900B CE	SpO <sub>2</sub> adapter cable (12-pin) Adapts 9-pin sensors to 12-pin sockets Length: 3m (9.8')	1 cable	Fetal/Maternal Monitor M1350B,M1350C; Defibrillators, M1722A, M1722B, M1723A, M1723B, M2475B; SpO <sub>2</sub> Module, M1020A,
2 M1940A CE	SpO <sub>2</sub> Adapter Cable (12-pin) Adapts 8-pin sensors to 12-pin sockets Length: 2m (6.6')	1 cable	
3 M1941A CE	SpO <sub>2</sub> extension cable (8-pin) Adds additional length to 8-pin sensors to 8-pin sockets Length: 2m (6.6')	1 cable	SpO <sub>2</sub> Module: M1020B, Opt. A01 Multi-measurement servers: M3000A, M3001A, M3002A, Opt. A01 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863069, 863071, 863072, 863073, 863074, 863077, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283
4 M1943A CE	SpO <sub>2</sub> adapter cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 1.1m (3.3')	1 cable	
4 M1943AL CE	SpO <sub>2</sub> adapter cable (8-pin) Adapts 9-pin sensors to 8-pin sockets Length: 3m (9.8')	1 cable	



# Pulse oximetry

## SpO<sub>2</sub> adapter cables



### SpO<sub>2</sub> adapter cables (continued)

Product #	Description	Packaging	Use with:
1 989803148221 CE	Masimo® MP10 adapter cable (LNCS sensors) Length: 3m (10')	1 cable	SpO <sub>2</sub> Modules: M1020B Opt. A01 and opt A03 Multi-measurement servers: M3001A, M3002A, Opt. A01 and Opt. A03 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863069, 863070, 863071, 863072, 863073, 863074, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283
2 M1020-61100 CE	Masimo® MP12 adapter cable (LNOP sensors) Length: 3.6m (11.8')	1 cable	SpO <sub>2</sub> Modules: M1020B Opt. A01 and opt A03 Multi-measurement servers: M3000A, M3001A, M3002A, Opt. A01 and Opt. A03 Patient Monitors: M8102A, M8105A, M8105AS, 863063, 863064, 863065, 863066, 863068, 863071, 863072, 863073, 863074, 863080, 863081, 863082, 863085, 863086, 863087, 863088, 863264, 863265, 863266, 863276, 863277, 863278, 863279, 863283
3 M1943NL CE	OxiMax® SpO <sub>2</sub> adapter cable (8-pin) Adapts 9-pin sensors to IntelliVue A02/A04 sockets Length: 3m (9.8')	1 cable	SpO <sub>2</sub> Modules: M1020B Opt. A02 and Opt. A04 Multi-measurement servers: M3001A, M3002A, Opt. A02 and Opt. A04 Patient Monitors: M8102A, M8105A, M8105AS, Opt. SP4



# Temperature



## Skin surface temperature probe

Product #	Description	Packaging	Use with:
1 21091A CE	Disposable sterile skin surface temperature probe Length: 30" (76cm) Shelf life: at least 12 months	1 box = 20 probes	Temperature module: M1029A; Multi-measurement servers: M3000A, M3001A, M3002A Patient monitors: 863064, 863065, 863066, 863068, M3015A, M8102A, M8105A, M3015A
2 21078A CE	Reusable, skin surface temperature probe, autoclaveable (2-pin connector) Diameter: 9.5 mm (0.37") Length: 3.1 m (10.2')	1 probe	Temperature module: M1029A; Multi-measurement servers: M3000A, M3001A, M3002A, M3015A Patient monitors: 863064, 863065, 863066, 863068, 860335, 860338, M8102A, M8105A
3 989803162641 CE	Reusable, skin surface temperature probe, autoclaveable (1/4-inch phone plug connector) Diameter: 9.5 mm (0.37") Length: 3.1 m (10.2')	1 probe	Series 400 temperature monitors



# Temperature

## Temperature adapter cables



### Temperature adapter cables for Philips monitors

Product #	Description	Packaging	Use with:
1 21082A CE	Reusable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 3m (10')	1 cable	Temperature module: M1029A; Multi-measurement servers: M3000A, M3001A, M3002A Patient monitors: 863064, 863065, 863066, 863068, M3015A, M8102A, M8105A, M3015A
1 21082B CE	Reusable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 1.5m (5')	1 cable	
2 989803163281 CE	Temperature probe adapter cable. Required to connect 989803162641 (1/4" phone plug) temperature probe to Philips and legacy HP/Agilent temperature monitors with 2-pin input connector. Length: 3m (10')	1 probe	Philips and legacy HP/Agilent monitors with 2-pin input connector

### Temperature adapter cables for series 400 compatible, non-Philips monitors

3 989803162601 CE	Reusable temperature probe adapter cable for Series 400 monitors Length: 1.5m (5')	1 cable	Series 400 temperature monitors with 1/4 inch phone socket
3 989803162591 CE	Reusable temperature probe adapter cable for Series 400 monitors Length: 3m (10')	1 cable	Series 400 temperature monitors with 1/4 inch phone socket
4 989803163271 CE	Temperature probe adapter cable. Required to connect 21078A (2-pin) temperature probe to Series 400 temperature monitors with 1/4" phone plug Length: 3m (10')	1 probe	Series 400 temperature monitors with 1/4 inch phone socket



## Jaundice management



- > **BiliChek non-invasive bilirubin assessment tool** **54**
- > **BiliCal calibration tips** **55**
- > **BilEclipse phototherapy patch** **55**
- > **WeeSpecs** **56**



# BiliChek non-invasive bilirubin assessment tool



Using light instead of a needle, BiliChek allows you to test bilirubin levels without a blood sample. Only the individual calibration tip touches the baby. Used for testing at the baby's forehead or sternum, this innovative handheld unit gives instant results that allow you to make immediate treatment decisions. BiliChek can be used before, during and after phototherapy, making it ideal for the newborn nursery and the NICU. Clinically proven, the tool virtually eliminates the costs and time associated with blood draws and laboratory services.

1

Product #	Description	Quantity per case	12NC #
1 1067554 €€	BiliChek system with wall power supply (North America only)	1	989805644871
1 1067561 €€	BiliChek system with desktop power supply (North America only)	1	989805644881
1 1037632 €€	BiliChek advanced system (Not available in the USA)	1	989805605031
1058220 €€	BiliChek 2 wall medical power supply	1	989805605931
1058221 €€	BiliChek 2 desktop medical power supply	1	989805605941
1062275	BiliChek manual (English, French, German, Italian, Spanish, Portuguese)	1	989805622461
1062276	BiliChek manual (Danish, Dutch, Greek, Swedish, Norwegian)	1	989805622441
1099142	BiliChek manual (Czech, Finnish, Kazakh, Traditional Chinese, Polish, Russian, Turkish, Portuguese, Brazilian Portuguese)	1	989805644411
1103448	BiliChek manual (Hungarian, Lithuanian, Romanian, Korean, Bahasa (Indonesian), Simplified Chinese)	1	989805646001

## Replacement parts

Product #	Description	Quantity per case	12NC #
1061601	BiliChek 2 lamp replacement kit	1	989805605991
1070439	BiliChek 2 protective sleeve	1	989805622341
1102103	BiliChek battery door replacement kit	1	989805645201
1002738 €€	BiliChek 2 battery	1	989805603471
1055523 €€	BiliChek 2 charge base	1	989805605901
1016618	BiliChek cap plug	5	989805603671



# BiliCal calibration tips / BilEclipse phototherapy patch



1

## BiliCal calibration tips

The BiliCal individual calibration tip allows for calibration of the BiliChek device prior to each measurement. Individual calibration tips are disposable plastic, reducing the risk of infection or possible cross contamination.



2

## BilEclipse phototherapy patch

The BilEclipse patch isolates a section of skin before, during and after phototherapy so transcutaneous bilirubin measurements can be taken with BiliChek.

### BiliCal calibration tips


Product #	Description	Quantity per case	12NC #
1 B800-50	BiliCal calibration tips	50	989805607841 €€

### BilEclipse phototherapy patch

Product #	Description	Quantity per case	12NC #
2 B800-11	BilEclipse phototherapy patch	50	989805607751 €€



# WeeSpecs



### WeeSpecs phototherapy mask

WeeSpecs' two-piece design is made from a soft, hypo-allergenic, light-blocking material that can be adjusted for a custom fit. The inside frame of the goggles seals around the skeletal socket reducing pressure and allowing the eyes to open and close naturally.

### WeeSpecs Supreme phototherapy mask

WeeSpecs Supreme incorporates an anatomically contoured goggle design and dark fabric that provide superb light-blocking protection without placing pressure on the eyes and is designed to fit even the smallest babies. WeeSpecs are held in place with a layer of breathable, slip-resistant foam that helps reduce movement.

### WeeSpecs phototherapy mask

Protection for baby's eyes, two-piece design

Product #	Description	Quantity per case	12NC #
1 99071-01	WeeSpecs preemie – fits head circumference <28cm (<11")	50	989805607631 €€
1 99071-02	WeeSpecs small – fits head circumference <28cm to 34cm (<11 - 13.4")	50	989805607641 €€
1 99071-03	WeeSpecs medium – fits head circumference >34cm (>13.4")	50	989805607651 €€

### WeeSpecs Supreme phototherapy mask

An easy-to-adjust strap offers a wide sizing range and a snug fit.

Product #	Description	Quantity per case	12NC #
2 1022894	WeeSpecs Supreme preemie – fits head circumference <26cm (<10.2")	50	989805604201 €€
2 1022481	WeeSpecs Supreme small fits – fits head circumference <26cm to 34cm (10.2" to 13.4")	50	989805604161 €€
2 1022896	WeeSpecs Supreme medium – fits head circumference >34cm (>13.4")	50	989805604211 €€





## Positioning



> <b>Bendy Bumper</b>	<b>58</b>	> <b>Squishon 2 and 3 / Squishon mattress</b>	<b>63</b>
> <b>SnuggleUp</b>	<b>59</b>	> <b>Gel-E Donut</b>	<b>64</b>
> <b>Frederick T. Frog</b>	<b>60</b>	> <b>WeePee diapers</b>	<b>65</b>
> <b>Prone Plus</b>	<b>61</b>	> <b>Pee-on absorbent underpad</b>	<b>65</b>
> <b>IncuCover</b>	<b>62</b>		
> <b>Infatherm Transport Mattress</b>	<b>62</b>		



# Bendy Bumper



These flexible positioning aids stay in place until reshaped – to promote containment and flexion, providing a reflex stimulus for extremity extension and recoil.

## Better Bendy (wipeable)

	Product #.	Description	Quantity per case	12NC #
2	1029244	B3 Better Bendy Bumper complete reusable unit with a permanent, wipeable urethane cover extra small: 71cm (28") long	6	989805604581
3	1029243	B3 Better Bendy Bumper complete reusable unit with a permanent, wipeable urethane cover small: 91.5 cm (36") long	6	989805604571

## Baby Bendy

	Product #.	Description	Quantity per case	12NC #
1	92065-B	Baby Bendy Bumper complete unit with rod designed for infants <1200g (<2.6lb)	6	989805606731
	92065-B1	Replacement foam body for Baby Bendy Bumper	6	989805606741

## Bendy Bumper

	Product #.	Description	Quantity per case	12NC #
1	92065-A	Bendy Bumper complete unit with rod designed for infants >1200g (>2.6lb)	6	989805606671
	92065-A1	Replacement foam body for Bendy Bumper	6	989805606681

## Covers

	Product #.	Description	Quantity per case	12NC #
	92065-B2	Replacement cloth covers for Baby Bendy Bumper or B3 extra small	12	989805606761
	92065-B3	Disposable covers for Baby Bendy Bumper or B3 extra small	50	989805606771
	92065-A2	Replacement cloth covers for Bendy Bumper or B3 small	12	989805606701
	92065-A3	Disposable covers for Bendy Bumper or B3 small	50	989805606711



# SnuggleUp



The SnuggleUp is a soft, cozy nest that helps to provide proper positioning and physiological stability for preterm and ill babies. Its padded footroll and soft, adjustable straps allow movement while gently maintaining appropriate positioning.

### Cloth style

SnuggleUps are made from a soft, washable cotton/polyester interlock knit that can be washed and reused.

Product #	Description	Quantity per case	12NC #
91015-X	Extra small for infants <900g (<2lb)	12	989805606311
91015-A	Small for infants 800 to 1350g (1.8 to 3lb)	12	989805606281
91015-B	Medium for infants 1300 to 2250g (2.8 to 5lb)	12	989805606291
91015-C	Large for infants > 2250g (>5lb)	12	989805606301

### Disposable style

Disposable SnuggleUps are particularly useful for admissions, procedures, transport, or when laundry facilities are not available.

Product #.	Description	Quantity per case	12NC #
97013-XDIS	Extra small for infants <900g (<2lb)	12	989805607421
97013-ADIS	Small for infants 800to 1350g (1.8 to 3lb)	12	989805607391
97013-BDIS	Medium for infants 1300 to 2250g (2.8 to 5lb)	12	989805607401
97013-CDIS	Large for infants > 2250g (>5lb)	12	989805607411



# Frederick T. Frog



1



2

Our versatile, frog-shaped beanbags, designed by a NICU Nurse, help achieve clinicians' goals for developmentally supportive positioning. Place Freddy around the head, neck, or hips. Or use it to position an extremity and provide tactile support. Available in two versions: washable/reusable, or wipeable/reusable with launderable or single-use cover.

## Frederick T. Frog – machine washable, dryable

Product #	Description	Quantity per case	12NC #
1 97040	Frederick T. Frog	6	989805603341

## Frederick T. Frog – wipeable

(use with launderable or disposable cover)

Product #	Description	Quantity per case	12NC #
2 1084677	Frederick T. Frog, wipeable with covers	6	989805635141
1084678	Launderable cover	50	989805635151
1084679	Disposable cover	50	989805635131



# Prone Plus



Provide ventral support to infants with a comfortable, hourglass-shaped pad. Prone Plus helps support newborns in the proper prone position by using gravity to facilitate natural rounding of the shoulders, hand-to-mouth coordination, and head-and-trunk alignment.

## Prone Plus

Provides standard and consistent prone support while helping to maintain a secure, balanced and flexed position.

Product #	Description	Quantity per case	12NC #
1034167	Prone Plus with knit cover extra small for infants <800g (<1.8lb)	6	989805604811
1034168	Prone Plus with knit cover small for infants 800 to 1250g (1.8 to 2.8lb)	6	989805604821
1034169	Prone Plus with knit cover medium for infants 1250 to 2000g (2.8 to 4.4lb)	6	989805604831
1034986	Prone Plus with knit cover large for infants >2000g (>4.4lb)	6	989805604891
1034170	Washable knit cover (extra small)	6	989805604841
1034171	Washable knit cover (small)	6	989805604851
1034172	Washable knit cover (medium)	6	989805604861
1034987	Washable knit cover (large)	6	989805604901



# IncuCover



The heart-print, heavy-duty fabric of the IncuCovers blocks out light and the neutral-colored underside is designed for accurate evaluation of skin color and perfusion. The full-length IncuCover has two long front flaps which can be drawn to create a dark environment.

IncuCovers provide babies with a darkened, protective environment.

Product#	Description	Quantity per case	12NC #
93029-04	Fits GE Giraffe® Omni-Bed®	3	989805606891

# Infatherm Transport Mattress



The Infatherm warming gel mattress is designed to gently support an infant's entire body. In the NICU, nursery, delivery room, ER or OR, it is used to provide relief at pressure points and assist in maintaining a neutral thermal environment.

## Infatherm Transport Mattress – disposable transport mattress with warming gel

When activated, the gel solution heats to 47°C/117°F, providing necessary thermoregulation. An infant can lie directly on its soft, non-slip fabric.

Product #	Description	Quantity per case	12NC #
1015	Disposable transport mattress	6	989805616831 €€



# Squishon 2 and 3 / Squishon Mattress



1



2



3

### Squishon 2 and 3

Provides an exceptionally soft, pliable surface to support and cradle a baby’s head and/or body. The rectangular shape allows for head movement while maintaining a supportive surface. Made from a soft material and filled with a non-toxic gel.

### Squishon Mattress

A gel-filled mattress provides a soft surface designed to help relieve pressure points, protect fragile skin and encourage proper positioning. Caregivers can reposition the baby and shift pressure points without having to handle the infant.

### Squishon 2 and 3 – versatile gel-filled cushions

Squishon 2 to support the extremely low birth weight infant.  
Squishon 3 for older preemies and babies with more resilient skin.

Product #	Description	Quantity per case	12NC #
1 91033-S	Squishon 2 15x22.5cm (6x9") <sup>1</sup>	12	989803196361 CE
1 91033-3	Squishon 3 15x22.5cm (6x9") <sup>1</sup>	12	989805606461 CE
1 91033-A	Disposable covers 18x23cm (7x9")	50	989805606481

### Squishon Mattress – pliable gel-filled mattress

Product #	Description	Quantity per case	12NC #
2 91017	Squishon mattress 30x55cm (12x22") <sup>1</sup>	1	989805603301 CE
2 91017-A	Disposable covers	25	989805606321

<sup>1</sup> Has an 18 month shelf life from the date of manufacture.



# Gel-E Donut



## Gel-E Donut

An exceptionally soft skin combined with a unique gel formulation make Gel-E Donuts soft, yet supportive. We designed it this way to help relieve pressure caused by prolonged immobility and other conditions where frequent handling may be contraindicated. Reengineered in 2015 to change to a non-porous skin.

Product #	Description	Quantity per case	12NC #
92025-XS	Extra small 18 cm diameter (7"), fill weight: 325g (11.5oz), purple <sup>1</sup>	12	989803196371 €€
92025-S	Small 18 cm diameter (7"), fill weight: 400g (14oz), green <sup>1</sup>	12	989803196381 €€
92025-M	Medium 18 cm diameter (7"), fill weight: 450g (16oz), peach <sup>1</sup>	12	989803196391 €€
91033-A	Disposable covers	50	989805606481

<sup>1</sup> Has an 18 month shelf life from the date of manufacture.





# WeePee diapers / Pee-on absorbent underpad



### WeePee diapers

WeePee diapers are specifically designed for premature infants to provide the right fit. The highly absorbent rectangular pad as well as the soft, flexible, non-adhesive closure system promotes developmentally appropriate positioning. Three sizes for a custom, comfortable fit in color-coded packaging to differentiate sizes at a glance.

### Pee-on absorbent underpad

The Pee-on is a 10" x 14" non-allergenic, non-toxic, appropriately-sized disposable underpad that helps to protect linens from urine, blood and spills. It can be used as a mattress protector or with all SnuggleUp positioning aids.

	Product #	Description	Quantity per case	12NC #
1	1017824	Extra small for infants <800g (<1.7lb)	360*	989805603721
1	1017855	Small for infants 800 to 1250g (1.7lb to 2.7lb)	360*	989805616941
1	1017856	Medium for infants 1250 to 2250g (2.7 to 5lb)	360*	989805603741
2	97047-B	Pee-on absorbent underpad	600	989805607451

\*30 diapers per sleeve; 12 sleeves per case



## Calming and soothing



> <b>Wee Thumbie</b>	<b>67</b>
> <b>Wee Soothie</b>	<b>67</b>
> <b>NICU Soothie / Soothie</b>	<b>68</b>
> <b>Super Soothie</b>	<b>69</b>
> <b>WubbaNub</b>	<b>69</b>
> <b>Sweet-Ease cup/vial</b>	<b>70</b>
> <b>Snoedel</b>	<b>71</b>



# Wee Thumbie / Wee Soothie



1



2



3

### Wee Thumbie

The Wee Thumbie is one of the first pacifiers designed specifically for use by very low birth weight babies. It is modeled after the size and shape of the preemie thumb, facilitating the important sucking behaviors normally learned in utero. The notched area allows space for CPAP tubing, mechanical ventilation, or tube feeding. The Wee Thumbie is not manufactured with DEHP or natural rubber latex.

### Wee Soothie

The Wee Soothie provides a superb fit for growing preemies or smaller term babies during their hospital stay. The notched area allows space for CPAP tubing, mechanical ventilation, or tube feeding. Our Wee Soothie pacifier is not manufactured with BPA, DEHP or natural rubber latex.

### Wee Thumbie

For babies <30 weeks gestation; for hospital use only

	Product #	Description	Quantity per case	12NC #
1	1046741	Wee Thumbie, aqua	100	989805605381
2	1046742	Wee Thumbie, purple	100	989805605391

### Wee Soothie

For babies 30–34 weeks gestation; for hospital use only

	Product #	Description	Quantity per case	12NC #
3	96003-N	Wee Soothie, natural	100	989805607121
3	96003-V	Wee Soothie, vanilla	100	989805607131



# NICU Soothie / Soothie



1

### NICU Soothie

The NICU Soothie is a premium pacifier designed for babies of 34 weeks gestation to 3 months old. The notched area allows for CPAP tubing, mechanical ventilation or an NG feeding. The NICU Soothie is not manufactured with BPA, DEHP or natural rubber latex.



2

### Soothie

The Soothie is a premium pacifier designed for newborns and babies up to 3 months old without teeth who are successfully breast or bottle feeding. The Soothie is not manufactured with BPA, DEHP or natural rubber latex.

### NICU Soothie

For babies >34 weeks gestation; for hospital or home use

	Product #	Description	Quantity per case	12NC #
1	1029232	NICU Soothie, natural (notched)	100	989805604551
1	1029233	NICU Soothie, vanilla (notched)	100	989805604561

### Soothie

For babies without teeth; for hospital or home use

	Product #	Description	Quantity per case	12NC #
2	96004-N	Soothie, natural	100	989805607151
2	96004-V	Soothie, vanilla	100	989805607191



# Super Soothie / WubbaNub



### Super Soothie

Similar to Soothie, but more durable, the Super Soothie is not manufactured with BPA, DEHP or natural rubber latex. Its one-piece construction adheres to AAP guidelines.

### WubbaNub

The WubbaNub uses the same Soothie pacifier that babies love and hospitals trust. This soft pacifier for newborns without teeth is specially designed to stay in place without clip or cord. A vibrantly colored, soft, plush animal gives the baby something to hold for better positioning. Available as a dog or duck, WubbaNub Soothies are not manufactured with BPA, DEHP, or natural rubber latex.

### Super Soothie

For babies with teeth; for hospital or home use

	Product #	Description	Quantity per case	12NC #
1	1031640	Super Soothie, pink (natural)	100	989805604751
1	1031639	Super Soothie, blue (natural)	100	989805604741

### WubbaNub

Pacifier with plush animal for babies without teeth; for hospital or home use. Available in the USA, Canada, Korea, Australia and New Zealand only.

	Product #	Description	Quantity per case	12NC #
2	99011-01	WubbaNub dog	1	989805607571
3	99011-02	WubbaNub duck	1	989805607581



# Sweet-Ease cup/vial

This 24% sucrose and purified water solution is a proven way to calm and soothe infants in the NICU, PICU, or well-baby nursery. It's available in a 2ml vial or 15ml cup to simplify putting it directly onto the tip of an infant's tongue.



1

### Sweet-Ease 15ml cup

This 15ml cup is aseptically packaged with a peel-off lid that is suitable for dipping a pacifier or for administration via a dropper or syringe.



2

### Sweet-Ease 2mL vial

This 2ml sterile vial includes 1D and 2D barcode stickers and simplifies clinician workflow with fewer steps and less clean-up.

## Sweet-Ease

	Product #	Description	Quantity per case	12NC #
1	99044	Sweet-Ease natural (15ml cup) <sup>1</sup>	200	989805603401
2	1093977	Sweet-Ease natural (2ml vial) <sup>1</sup>	96	989805637451

<sup>1</sup>Has a two-year shelf life from the date of manufacture



# Snoedel



The Snoedel, originating from a Dutch family tradition, offers parents a means to comfort and bond with their child even when the baby is in the hospital. The pure sheep wool incorporated into the high-quality flannel doll absorbs the parent's scent and slowly releases it, supporting the sensitive olfactory system of the infant.

## Snoedel

Snoedel flannel doll and bonding/positioning aid

Product #	Description	Quantity per case	12NC #
96029	Snoedel (6 blue/6 pink)	12	989805603331 €€
96029-B	Snoedel (12 unisex print)	12	989805607221 €€



## Therapy support products



- › **Nasal aspirators** **73**
- › **WeeSight transilluminator** **74**
- › **Infant Heel Warmer / Heel Snuggler** **75**
- › **Specialty feeding products by Pigeon** **76**
- › **Sali-Wipes** **77**
- › **Temperature Trend Indicators** **77**





# Nasal aspirators



1

## BBG nasal aspirator

Developed in collaboration with clinicians, the BBG nasal aspirator provides non-invasive suction of the nasal cavity for newborns, infants and children. Two sizes, for larger and smaller infants. A bulb-shaped interface provides an easy seal at the nare to protect the fragile nasal mucosa. BBG's versatility makes it ideal for the NICU, PICU, delivery room, emergency room and general pediatrics.



2

## Bulb-Style nasal aspirator

This aspirator provides gentle suction to remove excess mucus or secretions from a baby's nasal passage. The First Years® bulb-style aspirator, features an over-insertion guard and a pull-out plug that allows thorough cleaning. Philips is pleased to be the exclusive provider of this pediatrician-approved nasal aspirator to hospitals.

## Nasal aspirator

	Product #	Description	Quantity per case	12NC #
1	93032-S	BBG nasal aspirator 11.5cm (4.5") length	50	989805606961 CE
1	93032-02	BBG nasal aspirator 8.25cm (3.25") length	50	989805606951 CE
2	01053-1	Bulb-style nasal aspirator	24	989805606041



# WeeSight transilluminator



The WeeSight transilluminator helps to locate the tiny veins in a neonate's small limbs for improved insertion of an IV.



## WeeSight transilluminator

Uses two AA alkaline batteries (included). Output light wave length of approximately 629 nm.

Product #	Description	Quantity per case	12NC #
1017920	WeeSight transilluminator	1	989805603761 CE



# Infant Heel Warmer / Heel Snuggler



1



2

## Infant Heel Warmer

This easy-to-use heel warmer provides steady warming without skin adhesives or wet wash towels. The temperature of this heel warmer is dependable during each use providing consistency for the heel stick procedure.

## Heel Snuggler

The Heel Snuggler is the only infant heel warmer designed to conform truly to the infant's heel while staying firmly in place without skin adhesives. The unique design is shaped like a butterfly with two soft perforated straps for a comfortable and secure fit.

	Product #	Description	Quantity per case	12NC #
1	1223	Infant Heel Warmer	100 <sup>1</sup>	989805603201 CE
2	99047	Heel Snuggler	100 <sup>1</sup>	989805603411 CE

<sup>1</sup>25 per box; 4 boxes per case



## Specialty feeding products by Pigeon



The cleft palate nipple is specially designed to feed cleft palate infants and infants with poor sucking strength. The Y-cut nipple has a thick and thin sidewall with a one-way valve that prevents excessive air intake and allows milk to flow only when sucked by the baby. Not manufactured with latex. Available in the USA only.

The Pigeon food feeder bottle is a 120ml squeezable bottle with an attached spoon to assist with weaning or to use with babies who have difficulty feeding.

Product #	Description	Quantity per case	12NC #
98072	Food feeder bottle	1	989805603351
1104725	Cleft palate bottle with 2 nipples	10	989805646671
1104727	Cleft palate nipple, regular size	5	989805646691
1104758	Cleft palate nipple, small size	5	989805646701



# Sali-Wipes / Temperature Trend Indicators



1



2



## Sali-Wipes

Sali-Wipes sterile saline (isotonic) wipes come in convenient foil packs. Nursing organizations recommend use of sterile isotonic saline for cleansing, as alcohol wipes can be damaging to fragile infant skin and may be absorbed into the bloodstream. Available in the USA only.

## Temperature Trend Indicators

Temperature Trend Indicators provide rapid, accurate temperature information non-invasively for consistent thermal monitoring of patients. Available in the USA only.

Product #.	Description	Quantity per case	12NC #
1 98098-013	Sali-Wipes 10cm x 15.2cm (4" x 6")	576	989805607531

## Temperature Trend Indicators

Product #.	Description	Quantity per case	12NC #
2 LC105D	LC Accu-Bar Temperature Trend Indicator and dispenser box (one degree thermometer) 94F-105F	100	989805623001



## Education and consulting programs



› **Education and consulting programs**

**79**



### NICU Core Development Program

Nurses who support the neuroprotective care and specialized needs of fragile newborns in the NICU require special training to work safely and effectively. The NICU Nurse Core Development Program provides nurses with the knowledge, skills, and competencies needed to deliver responsive care to neonatal patients and their families. This comprehensive program is customized to support an institution's educational objectives, clinical roles, job tasks, workflows, policies and procedures. It can be used to onboard new neonatal nurses, and it provides a learning environment for on-going professional development for staff.

### Wee Care Neuroprotective NICU Program

Wee Care is a customized, internationally recognized education and quality improvement program designed to enhance your goals for family-centered neuroprotective Developmental Care practices in the NICU and beyond. It combines evidence-based best practices around seven core measures for infant development. Through change management and education, Wee Care helps to promote natural development, healthy growth, and smooth transitions out of the NICU to the home.

### Preemie for a Day

This popular, interactive, multi-sensory training program contrasts the experience of an infant in a traditional NICU environment to that of an infant in a developmentally supportive NICU. The practice areas addressed are admissions, positioning and handling, feeding, and family and staff partnerships.

### Neonatal Integrative Developmental Care Model



© Koninklijke Philips N.V., 2016. All rights reserved.



### **Cue-Based Feeding: Implications for Clinical Practice**

This innovative program is aimed at individualizing the feeding experience for infants in the NICU through the practice of infant-driven feeding. The program is designed to address the physiologic, psychological, developmental and environmental issues that impact feeding experiences for the preterm and/or critically ill- hospitalized infant. You'll gain a deeper understanding of the crucial nature of the feeding experience and the implications of feeding interventions that impact a vulnerable infant. You'll come away with tools to promote and implement change in feeding practices within your own unit.

### **Wee Care Sacred Hour Maternity Program**

Wee Care Sacred hour maternity program is an educational and transformational change

consulting engagement that supports the hospital's efforts in neuroprotective family-centered Developmental Care in the perinatal areas (L&D/PPM/NB/NICU). The goal is to promote skin-to-skin contact (SSC) during the first hour after birth to increase early breastfeeding initiation and exclusive breastfeeding rates.

### **Wee Focus**

The Wee Focus program draws from the success of the internationally recognized Global Wee Care® Program to "deep-dive" into a single core measure from the evidence-based Philips Neonatal Integrative Developmental Care model. Wee Focus is a process improvement program designed to help individual teams implement evidence-based best practices aimed at achieving a specific initiative. This 10-14 week clinical transformation performance improvement program is tailored to the unit's specific

goals and improvement needs based on a single core measure. The Wee focus program is aimed at implementing neuroprotective family-centered Developmental Care practices in the NICU. The goal is to enhance the NICU experience for newborns and their families through education and change management.

### **Late Preterm Infant**

The Late-Preterm Infant (LPI) program shares insights relating to the complexity of the late preterm infant and provide caregivers with a better understanding of the hidden challenges these infants face. Program topics include recognizing physiologic and behavioral differences between the preterm, late preterm, and full-term infant, identifying the LPI at risk for complications, and fostering developmentally supportive care throughout admission and the hospital stay.

Learn more about these and new Mother and Child Care education programs at:

<http://www.learningconnection.philips.com/en/catalog/product-group/mother-and-child-care>



